

ASMS 2014, (Jun 16, 2014, Baltimore, USA)

# Development of a novel amino acids analysis column for LC-MS without derivatization

YAZWA Itaru (yazawa@imtakt.com) and TACHIKAWA Hiroshi

Imtakt Corp.

Kyoto Research Park, Kyoto 600-8813, Japan



# OVERVIEW

We have developed a novel amino acid separation column for LC-MS without derivatization.

- Mixed-mode stationary phase with NP+IEX
- Achieved to separate and detect or non-label 55 amino acids
- Simple gradient separation with two mobile phases
- High-throughput separation within 1 to 10 min
- Successful separation for isomers (Leu/Ile, GABA etc.)
- Good reproducibility and linearity
- Applicable to dipeptides analysis
- Extensibility to various amino acid related compounds

This novel HPLC method will be a powerful tool for amino acid LC-MS analysis in many different biochemistry applications.

# INTRODUCTION

There are four established methods for analyzing amino acids: pre-labeled, post-labeled, ion-pairing reversed-phase, and normal-phase, but each of these methods has disadvantages.

The pre-labeled method has problems with derivitization efficiency and cost, while the post-labeled method is usually not compatible with LC-MS due to non-volatile mobile phases. The ion-pairing reversed-phase method has difficulty separating polar amino acids; on the other hand, the normal-phase mode has problems separating all the compounds, especially the Leu and Ile isomers.

We have developed a novel amino acid separation column for LC-MS(/MS) which can separate all 20 amino acids in protein using a mixed-mode stationary phase structure. We have also estimated separation and detection characteristics using LC-MS instruments.

## METHOD

We used a novel amino acid analysis column for LC-MS. This stationary phase consists of 3µm silica particles modified with ion ligands which provide IEX and normal phase mixed-mode separation.

LC-MS instrument:

A mass spectrometer, LCMS-2020 (Single Quad) with LC-30A binary pumps (Shimadzu Corporation, Kyoto)

Reagents:

All of reagents including standard amino acids and bovine serum were purchased as commercial reagents.

## RESULTS AND DISCUSSION

Described in each figure.


## CONCLUSION

This novel HPLC method will be a powerful tool for amino acid LC-MS analysis in many different biochemistry applications.



# LC-MS analysis of 20 amino acids

## Single Quad. MS




20 standard amino acids of protein were separated and detected on LC-MS (Single Quad).  
 Leucine and Isoleucine isomers were also separated in several minutes.



## LC-MS/MS (MRM) Analysis

**Triple Quad. MS**  
**100 nmol/mL**




20 standard amino acids of protein were separated and detected on triple quad. (MRM) with excellent sensitivity.



## Simple Gradient Conditions

Single Quad. MS  
100 nmol/mL



Simple gradient elution condition was also successful for standard amino acids separation in 10min with low pressure.



## 5min High Throughput Analysis

**Amino acids standard mixture**  
(100 nmol / mL)

**Intrada Amino Acid, 150 x 2 mm**


**A: tetrahydrofuran / 20 mM ammonium formate / formic acid**  
= 75 / 25 / 0.3

**B: 200 mM ammonium formate / acetonitrile = 80 / 20**

0 %B (0-2.5 min), 100 %B (2.5-5 min)

0.6 mL/min (45 MPa), 37 deg.C, 5 uL

ESI (SIM, positive)



High throughput analysis using UHPLC system was achieved to analyze 20 amino acids within 5min.



5



## Accurate method with 100 x 3 mm

### 100 x 3 mm

*Intrada Amino Acid, 100 x 3 mm*

*A: ACN /THF /25 mM HCOONH<sub>4</sub> /HCOOH = 9 /75 /16 /0.3*

*B: ACN /100mM HCOONH<sub>4</sub> = 20 /80*

*0 %B (0-3 min)*


*0-17 %B (3-9 min)*

*100 %B (9-16 min)*

*0.6 mL/min (9MPa)*

*35 deg.C, 5uL (0.1N HCl)*

*ESI (SIM, positive)*



**100 x 3 mm column provided more accurate analysis for 20 standard amino acids with better peak shape and separation.**

6



# One-Minute High Throughput Analysis

*Intrada Amino Acid, 10 x 2 mm*

*A: acetonitrile / formic acid = 100 / 0.1*


*B: 100 mM ammonium formate*

*15-100 %B (0-0.8 min), 100 %B (0.8-1.0 min)*

*0.4 mL/min (1.6 MPa), 35deg.C, 1 uL (0.1N HCl)*

*ESI (SIM, positive)*


**10 x 2 mm**



**Amino acid related compounds were analyzed by one-minute high throughput separation with 10mm length small column.**

**7**

## Batch-to-Batch Reproducibility (20 AA)

**150 x 3 mm**

**Intrada Amino Acid, 150 x 3 mm**





**A: ACN / THF / 25mM HCOONH<sub>4</sub> / HCOOH = 9 / 75 / 16 / 0.3 (v/v/v/v)**

**B: ACN / 100mM HCOONH<sub>4</sub> = 20 / 80 (v/v)**

**0 %B (0-12 min), 0-17 %B (12-26 min), 100 %B (26-45 min)**

**0.4 mL/min (7MPa), 35 deg.C, 1.2 uL (1 umol/mL in 0.1N HCl), ESI, positive**

**Amino acid analysis packing materials were showed excellent batch-to-batch reproducibility.**



*Intrada Amino Acid, 50 x 3 mm, LCMS-2020 (single quad.), ESI (SIM, positive)*

**Typical standard (aromatic, aliphatic, acidic and basic) amino acids showed good results for linearity and sensitivity.**



## Repeatability

**Amino Acid**    **%RSD(t<sub>R</sub>)**    **%RSD(Area)**

<b>Trp</b>	<b>0.48</b>	<b>0.82</b>
<b>Phe</b>	<b>0.49</b>	<b>0.97</b>
<b>Tyr</b>	<b>0.70</b>	<b>1.09</b>
<b>Met</b>	<b>0.45</b>	<b>1.21</b>
<b>Leu</b>	<b>0.43</b>	<b>2.33</b>
<b>Ile</b>	<b>0.39</b>	<b>3.11</b>
<b>Val</b>	<b>0.41</b>	<b>2.26</b>
<b>Glu</b>	<b>0.40</b>	<b>1.71</b>
<b>Pro</b>	<b>0.34</b>	<b>2.08</b>
<b>Asp</b>	<b>0.28</b>	<b>2.60</b>
<b>Thr</b>	<b>0.38</b>	<b>2.39</b>
<b>Ala</b>	<b>0.24</b>	<b>3.10</b>
<b>Ser</b>	<b>0.22</b>	<b>2.71</b>
<b>Gln</b>	<b>0.22</b>	<b>3.21</b>
<b>Gly</b>	<b>0.21</b>	<b>2.92</b>
<b>Asn</b>	<b>0.20</b>	<b>2.35</b>
<b>(Cys)2</b>	<b>0.04</b>	<b>4.14</b>
<b>His</b>	<b>0.08</b>	<b>4.30</b>
<b>Lys</b>	<b>0.08</b>	<b>4.85</b>
<b>Arg</b>	<b>0.11</b>	<b>4.11</b>

*Intrada Amino Acid, 50 x 3 mm*

*A: ACN /THF /25mM HCOONH<sub>4</sub> /HCOOH = 9 /75 /16 /0.3*

*B: ACN /100 mM HCOONH<sub>4</sub> = 80 / 20*

*0 %B (0-2.5 min)*

*0-17 %B (2.5-6.5 min)*

*100 %B (6.5-10 min)*

*0.6 mL/min*

*35 deg.C*

*5uL (500nmol/mL)*

*ESI, Positive*

*n = 15*

*ISTD: theanine*

**Amino acid analysis column showed excellent repeatability in retention and acceptable for peak area.**

Intrada Amino Acid, 50 x 3 mm

A: acetonitrile /tetrahydrofuran /25mM ammonium formate /formic acid  
= 9 / 75 / 16 / 0.3

B: acetonitrile / 100mM ammonium formate = 20 / 80


0 %B (0-2.5min)

0-17 %B (2.5-6.5 min)


100 %B (6.5-10 min)

0.6 mL/min (6MPa), 35 deg.C, 5 uL (1000nmol/mL)

LCMS-2020 (Single quad.), ESI (SIM, positive)




55 standard amino acids were separated and detected on LC-MS. Almost amino acids were analyzed within 10 minutes.




***Leucine isomers, same m/z value, were well separated in short time.***


*m/z 104.1*



**3-aminobutyric acid**  
**(beta-aminobutyric acid)**




**4-aminobutyric acid**  
**(GABA, gamma-aminobutyric acid)**




**2-aminobutyric acid**  
**(alpha-aminobutyric acid)**

*Intrada Amino Acid, 75 x 2 mm*  
*A: acetonitrile / formic acid = 100 / 0.3*  
*B: 100 mM ammonium formate*  
*15-100 %B (0-8 min), 15 %B (8-10 min)*  
*0.3 mL/min (4 MPa), 35 deg.C*  
*5 uL (10 ug/mL, 0.1N-HCl aq.)*  
*ESI (SIM, positive, m/z 104.1)*



**Not only alpha-amino acids but also beta- and gamma-amino acids isomers were well separated.**



**L-cystine ((Cys)2)**

**Intrada Amino Acid, 75 x 2 mm**

**A: acetonitrile / formic acid = 100 / 0.3**






**B: 100 mM ammonium formate**

**20-100 %B (0-10 min), 100 %B(10-12 min)**

**0.3 mL/min (5 MPa), 35 deg.C**

**2 uL (0.02-0.32 ug, 0.1N-HCl aq.)**

**ESI (SIM, positive)**

**N-(L-alpha-glutamyl)  
-L-glutamic acid  
(Glu-Glu)****L-anserine****L-carnosine****N-L-histidyl-L-leucine  
(His-Leu)**

**The novel amino acid analysis column could be applied for dipeptides.**



# Amino acids in human serum (HClO<sub>4</sub>)


Human Serum (commercial reagent)

← 0.4N HClO<sub>4</sub> (1:1, vol)  
← mixing, centrifugation

Supernatant

← 0.2µm filtration

Filtrate (Injection Sample)



The novel amino acid analysis column was applied for amino acids in serum and almost amino acids were detected.



## Determination of Neuronal Peptides and Catecholamines and Effects of Diethylstilbestrol on Male Rat Pituitary

Naoyuki Maeda<sup>1,2</sup>, Emi Tanaka<sup>1</sup>, Kanae Masu<sup>1</sup>, Kanako Okumura<sup>2</sup>, Yuki Ikeda<sup>2</sup>, Taku Miyasho<sup>2</sup>, Satoko Haeno<sup>2</sup> and Hiroshi Yokota<sup>2</sup>.

1) Japan Meat Science & Technology Institute

2) Rakuno Gakuen University

ACMS2013 (61st), Sept 10, 2013, Tsukuba, Japan

The separation was achieved using an Intradra Amino Acid  
(100 x 3 mm 3µm particle size, Imtakt)

(A) CH<sub>3</sub>CN / THF / 25mM HCOONH<sub>4</sub> / HCOOH = 10 / 80 / 10 / 0.4

(B) CH<sub>3</sub>CN/100mM HCOONH<sub>4</sub> = 20/80

0 %B (0-1 min)



0-17 %B (1-6.5 min)

100 %B (6.5-10 min)

100 %B hold (10-18 min)

400 µL/min flow rate at 35°C

5 µL injection



**Determination procedure for neuronal peptides and catecholamines were developed by MS analysis.**

**The regulatory factors for prolactin induction mediated by diethylstilbestrol were analyzed by LC-MS analysis using the selected reaction monitoring (SRM). Dopamine suppressing prolactin secretion did not decreased, on the other hand, vasoactive intestinal peptide (VIP) which mediates the acute release of prolactin was increased in the pituitary glands of the DES- treated rats.**