


Presto FF-C18

50 x 4.6 mm
150 x 4.6 mm


Application

Lipopolysaccharide (from *E.coli* O127)

リポ多糖(大腸菌O127由来)



Presto FF-C18, 50 x 4.6 mm
A: water / formic acid = 100 / 0.1
B: acetonitrile / formic acid = 100 / 0.1
0-70 %B (0-5 min)
0.5 mL/min (7 MPa), 37 deg.C
ELSD (spray chamber 30 deg.C,
drift tube 80 deg.C)
5 uL (10 ug)



Presto FF-C18, 150 x 4.6 mm
A: water / formic acid = 100 / 0.1
B: acetonitrile / formic acid = 100 / 0.1
0-70 %B (0-25 min)
0.2 mL/min (7 MPa), 37 deg.C
ELSD (spray chamber 30 deg.C,
drift tube 80 deg.C)
5 uL (10 ug)