


Unison UK-C1  
 Unison UK-Phenyl Cadenza CX-C18  
 Unison UK-C8 Cadenza CW-C18  
 Unison UK-C18 Cadenza CD-C18

100 x 4.6 mm


HPLC COLUMN APPLICATION

## 2,4,6-Tri ([1,1'-biphenyl] -3-yl) -1,3,5-triazine (Organic light-emitting diode material)

2,4,6-トリ ([1,1'-ビフェニル]-3-イル) -1,3,5-トリアジン  
 (有機エレクトロルミネッセンス材料)



2,4,6-tri ([1,1'-biphenyl] -3-yl) -1,3,5-triazine



100 x 4.6 mm

A: water / acetonitrile / formic acid = 50 / 50 / 0.1

B: acetonitrile / formic acid = 100 / 0.1

0-100%B (0-10min), 100%B (10-30min)

1.0mL/min (8MPa), 37°C, 260nm

10µL (25µg/mL, THF / formic acid = 100 / 1)