



Unison UK-Silica
 Unison UK-C8
 Unison UK-C18
 Cadenza CL-C18
 Cadenza CD-C18


150 x 3 mm
 75 x 3 mm

Application

Coenzyme Q10
 コエンザイムQ10



Unison UK-Silica
 150 x 3 mm
 hexane / IPA = 998 / 2
 0.5 mL/min, 37 deg.C
 270 nm, 2ug (2uL)



75 x 3 mm
 0.5 mL/min, 37 deg.C
 270 nm, 2ug (2uL)