


Scherzo SW-C18
 Scherzo SM-C18
 Scherzo SS-C18

100 x 4.6 mm


HPLC COLUMN APPLICATION

Tris (2,2'-bipyridyl) ruthenium (II) (Dye for dye sensitized solar cell)

トリス(2,2'-ビピリジル)ルテニウム(II) (色素増感太陽電池用色素)



tris (2,2'-bipyridyl) ruthenium (II)
 chloride hexahydrate



100 x 4.6 mm
 1.0 mL/min (9MPa), 37°C, 260nm